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## TEM staining report for cells and bacteria

### 1 Apparatus and reagents

#### 1.1 Major apparatus

Name	Producer	Model
Ultra microtome	Leica	Leica UC7
Diamond slicer	Daitome	Ultra 45°
Transmission Electron Microscope	HITACHI	HT7800
2000 meshes cuprum grids with formvar film	emcn	BZ11262a

#### 1.2 Major reagents

Name	Producer	Code
Fixative for TEM	bioosci	BP0130
Ethanol	Sinaopharm Group Chemical Reagent Co. LTD	100092183
Aceton	Sinaopharm Group Chemical Reagent Co. LTD	10000418
EMBed 812	SPI	02660-AB
OsO4	Ted Pella Inc	02602-AB
uranyl acetate	Ted Pella Inc	02624-AB
Agarose	Solarbio	A8201
Lead Citrate	emsdiasum	17800

### 2 Procedure

**2.1 Harvest samples and fixation:** Collect cells or bacteria precipitation after centrifuge. The TEM fixative was added to the tube and let the precipitation re-suspended in the fixative, and then fixed at 4°C for preservation and transportation.

**2.2 Agarose pre-embedding:** The fixed cells and bacteria were centrifuged. The 0.1 M PB (pH 7.4) was added into the tube after supernatant was discarded, and then the precipitation was re-suspended and washed in PB for 3min. This washing step was repeated for 3 times. The 1% agarose solution was prepared by heating and dissolving in advance. After being cooled, the agarose solution was added into the EP tube. Before agarose

solidification, the precipitation was suspended with toothpick and wrapped in the agarose. The size of cells precipitation should be no more than 1 mm<sup>3</sup>

**2.3 Post-fix:** Agarose blocks with samples avoid light post fixed with 1% OsO<sub>4</sub> in 0.1 M PB (pH 7.4) for 2 h at room temperature. After remove OsO<sub>4</sub>, the tissues are rinsed in 0.1 M PB (pH 7.4) for 3 times, 15 min each.

**2.1 Dehydrate at room temperature as followed:**

30% ethanol for 10 min;  
50% ethanol for 10min;  
70% ethanol for 10min;  
80% ethanol for 10min;  
95% ethanol for 10min;  
Two changes of 100% ethanol for 10 min;  
Finally two changes of acetone for 10min.

**2.4 Resin penetration and embedding as followed:**

Acetone: EMBED 812= 1:1 for 2-4h at 37°C

Acetone: EMBED 812= 1:2 overnight at 37°C

pure EMBED 812 for 5-8 h at 37°C;

Pour the pure EMBED 812 into the embedding models and insert the tissues into the pure EMBED 812, and then keep in 37°C oven overnight.

**2.5 Polymerization:** The embedding models with resin and samples were moved into 60°C oven to polymerize for more than 48h. And then the resin blocks were taken out from the embedding models for standby application at room temperature.

**2.6 Ultrathin section:** The resin blocks were cut to 60-80nm thin on the ultra microtome, and the tissues were fished out onto the 200 meshes cuprum grids with formvar film.

**2.7 Staining:** 2% uranium acetate saturated alcohol solution avoid light staining for 10 min, then rinsed in ultra pure water for 3 times. Lead citrate avoid CO<sub>2</sub> staining for 10 min, and then rinsed with ultra pure water for 3 times. After dried by the filter paper, the cuprum grids were put into the grids board and dried overnight at room temperature.

**2.8 Observation and images capture:** The cuprum grids are observed under TEM and take images.