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Experimental Report on Decalcification of Bone Tissue (EDTA)

1. Instruments and key reagents

1.1 Instruments

Instrument	Manufacture	Specifications/Model
Tissue processor	Changzhou Zhongwei Electronics Co., Ltd	TSJ-SD
Tissue embedder	Changzhou Zhongwei Electronics Co., Ltd	BMJ-A
Embedding frame	Fujian Qisheng Experimental Equipment Technology Co., Ltd	58444
Microtome	ThermoFisher Scientific	SHANDON FINESSE 325
Paraffin	Guangdong Dachuan Special Wax Co., Ltd	
Freezing table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Water bath - Slide drier	Changzhou Zhongwei Electronics Co., Ltd	PHY-III
Slide	Hubei BIOSSCI Biotech Co., Ltd	BP0510
Mixer	DLAB Scientific Co., Ltd	OS20-Pro
Analytical balance	Shanghai Yaoxin Electronic Technology Co., Ltd	LQ-A20002

1.2 Key reagents

Reagent	Manufacture	Specifications/Model
EDTA-2Na	Sinopharm	10009717
NaOH	Sinopharm	10019718

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2. Procedures

2.1 The fresh tissue was fixed with fix solution for more than 24 hours. The fix solution was poured out of the EP tube containing tissue after fixation. Then, the EP tube was filled with the prepared decalcifying solution and sealed.

2.2 Decalcification: The decalcified tissue shall be labeled and the temperature shall be generally controlled between 25°C and 30°C. The decalcified solution shall be replaced every 3 days. A small needle can be used to stab bone tissue (not target area), which can be easily pierced, indicating that decalcification is complete.

2.3 Tissue processing: The decalcified bone tissue was taken out from the tube and put into the embedding frame with corresponding label. Tissues were washed with running water and then put into tissue processor for dehydration and wax immersion.

2.4 Paraffin embedding: The wax-soaked tissue was embedded by using tissue embedder. Firstly, melted paraffin was put in the bottom of the embedded mold. Then, tissue was taken out of the embedding frame when the paraffin solidified slightly and put into embedded mold according to the requirements of the different section. Tweezers were used to gently press the tissue to make the tissue completely flat at the bottom of mold. Finally, the mold with tissue in it was covered by embedding frame with corresponding identifier and moved to the -20°C freezing table for cooling. The paraffin block was taken out from the mold after paraffin solidified.

2.5 The trimmed paraffin block was put on the microtome to be sliced with 4μm thickness. Sections were floated on warm water (42°C) to be flatten and were picked up vertically by using a glass slide. After tissue was dried, sections were put into incubator at 60°C for 30 minutes to 1 hour. Sections were taken out of the incubator and stored at room temperature.

3. Attention

The target area cannot be punctured. The needle head should be small. If the bone tissue has skin and muscle adhesion, it can be stripped before decalcification.

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4. Technical support

The technical support of decalcification experiment is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

Tel: 400 118 0100

Fax: +86-027-87382710

E-mail: support@biossci.com

Website: www.biossci.com