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## Experimental Report on Masson Staining

### 1. Instruments and key reagents

#### 1.1 Instruments

Instrument	Manufacture	Specifications/Model
Tissue processor	Changzhou Zhongwei Electronics Co., Ltd	TSJ-SD
Tissue embedder	Changzhou Zhongwei Electronics Co., Ltd	BMJ-A
Microtome	ThermoFisher Scientific	SHANDON FINESSE 325
Freezing table	Wuhan Junjie Electronics Co., Ltd	JB-L5
Water bath - Slide drier	Changzhou Zhongwei Electronics Co., Ltd	PHY-III
Slide	Hubei BIOSSCI Biotech Co., Ltd	BP0510
Upright microscope	Olympus	CX-31
Digital scanner	HAMAMATSU PHOTONICS	NanoZoomer®S360

#### 1.2 Key reagents

Reagent	Manufacture	Specifications/Model
Ethanol	Sinopharm	100092683
Clearer	Wuhan Tongsheng Technology Development Co., Ltd	
Neutral balsam	Wuhan Tongsheng Technology Development Co., Ltd	
Bouin's fix solution	Hubei BIOSSCI Biotech Co., Ltd	BP0140
Masson staining kit	Hubei BIOSSCI Biotech Co., Ltd	BP0310

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---

2.1 Deparaffinization and rehydration.

2.1.1 Tissue sections were immersed in clearer for 10min. Repeat this step two times, gently shaking off excess liquid between each step.

2.1.2 Tissue sections were immersed in progressively more dilute ethanol solutions and ultimately immersed in distilled water to rehydrate the tissue: Absolute ethanol for 5min, 95% ethanol for 5min, 85% ethanol for 5min, 75% ethanol for 5min. Rinsing with distilled water for 1min.

2.1.3 The dehydrated tissue sections were immersed in Bouin's solution or Zenker's solution overnight, then, rinsed with running water.

2.2 Sections were stained with hematoxylin solution (Harris) or iron hematoxylin for 5-10 minutes and slightly washed with running water.

2.3 Sections were differentiated with 0.8% - 1% hydrochloric acid alcohol and washed with running water for several minutes. Sections can also be treated with lithium carbonate solution to be bluer and washed with running water.

2.4 Sections were stained with ponceau acid fuchsin solution for 5-10 minutes and washed with running water.

2.5 Sections were treated with phosphomolybdic acid solution for about 5 minutes and then stained with aniline blue solution for 5 minutes without washing.

2.6 Sections were treated with 1% glacial acetic acid for 1 minute and dehydrated with 95% alcohol for several times.

2.7 The tissue sections were dehydrated with absolute alcohol and transparent with xylene, then mounted with neutral balsam. Examination with microscope.

### 3. Results

Collagen fibers were blue (counterstained with aniline blue solution) or green (counterstained with bright green solution). The cytoplasm, muscle fibers and red blood cells are red. Nuclei are blue brown.

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---

4.1 It is better to fix the tissue with Bouin's solution or Zenker's solution. If it has been fixed with 10% formaldehyde solution, the tissue sections can be dewaxed to water, and then put into Bouin's solution overnight at room temperature or put into a 37°C incubator for 1-2 hours. Then, washed sections with running water until the yellow disappeared before staining.

4.2 It need to be controlled under microscope during treatment of phosphomolybdic acid. It can be seen that muscle fibers are red and collagen fibers are light red.

4.3 It is important to control the time of differentiation with glacial acetic acid, if the differentiation is excessive, the blue of collagen will be too light. If the differentiation is insufficient, it is easy to overlap with red and become purple blue.

### 5. Technical support

The technical support of Masson staining experiment is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

Tel: 400 118 0100

Fax: +86-027-87382710

E-mail: support@biossci.com

Website: [www.biossci.com](http://www.biossci.com)