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## Experimental Report on Paraffin embedding

### 1. Instruments and key reagents

#### 1.1 Instruments

| Instrument               | Manufacture   | Specifications/Model |
|--------------------------|---|----------------------|
| Tissue processor         | Changzhou Zhongwei Electronics Co., Ltd                   | TSJ-SD               |
| Tissue embedder          | Changzhou Zhongwei Electronics Co., Ltd                   | BMJ-A                |
| Embedding frame          | Fujian Qisheng Experimental Equipment Technology Co., Ltd | 58444                |
| Microtome                | ThermoFisher Scientific                                   | SHANDON FINESSE 325  |
| Paraffin                 | Guangdong Dachuan Special Wax Co., Ltd                    |                      |
| Freezing table           | Wuhan Junjie Electronics Co., Ltd                         | JB-L5                |
| Water bath - Slide drier | Changzhou Zhongwei Electronics Co., Ltd                   | PHY-III              |
| Slide                    | Hubei BIOSSCI Biotech Co., Ltd                            | BP0510               |

#### 1.2 Key reagents

| Reagent | Manufacture                                     | Specifications/Model |
|---------|---|----------------------|
| Ethanol | Sinopharm                                       | 100092683            |
| Clearer | Wuhan Tongsheng Technology Development Co., Ltd |                      |

### 2. Procedures

2.1 The fresh tissue was fixed with fix solution for more than 24 hours. The tissue was removed from the fix solution after fixation. Then, the tissue was smoothed by

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---

using a scalpel in the fume hood and put into the embedding frame.

2.2 Washing: The tissue was put into tap water and washed for 20 minutes.

2.3 Tissue processing: Embedding frames containing tissue were put into the cassette and dehydrated with gradient ethanol in turn. In 75% ethanol for 2h, 85% ethanol for 2h, 95% ethanol(1) for 1h, 95% ethanol(2) for 2h, absolute ethanol(1) for 1h, absolute ethanol(2) for 1h, xylene(1) 1h, xylene(2) 1h, melted paraffin(1) 1h, melted paraffin(2) 1h.

2.4 Paraffin embedding: The wax-soaked tissue was embedded by using tissue embedder. Firstly, melted paraffin was put in the bottom of the embedded mold. Then, tissue was taken out of the embedding frame when the paraffin solidified slightly and put into embedded mold according to the requirements of the different section. Tweezers were used to gently press the tissue to make the tissue completely flat at the bottom of mold. Finally, the mold with tissue in it was covered by embedding frame with corresponding identifier and moved to the -20°C freezing table for cooling. The paraffin block was taken out from the mold after paraffin solidified.

2.5 The trimmed paraffin block was put on the microtome to be sliced with 4μm thickness. Sections were floated on warm water (42°C) to be flatten and were picked up vertically by using a glass slide. After tissue was dried, sections were put into incubator at 60°C for 30 minutes to 1 hour. Sections were taken out of the incubator and stored at room temperature.

### 3. Technical support

The technical support of paraffin section production is provided by Hubei BIOSSCI Biotech Co., Ltd (Wuhan Changyan Pathology technology Co., Ltd).

Tel: 400 118 0100

Fax: +86-027-87382710

E-mail: support@biossci.com

Website: www.biossci.com